



Eye Safety Filters

Features and Benefits

- 532 nm and 810 nm Eye Safety Filters for Endophotocoagulation
- Superior color balance and tissue visualization with 100 coating layers
- Compatible with standard Zeiss, Leica, Topcon, Wild, Moller Wedel microscopes
- 2-position, 4-channel fixed, 2-channel models available
- V notch filter design, low profile housing



Specifications are subject to change without notice. IRIDEX and the IRIDEX logo are registered trademarks and TruFocus is a trademark of IRIDEX Corporation.

Products are covered by one or more of the following U.S. patents: 5,085,492; 5,088,803; 5,372,595; 5,511,085; 5,982,789; 6,327,291; 6,540,391; and 6,733,940.

Other U.S. and International patents pending.

Clinical references available upon request.



IRIDEX

Our name is our reputation™