




ENT Laser Delivery Devices

Increasing Your Clinical Options

- For delivery of laser energy when performing ENT procedures
- Indicated for, but not limited to, incision, excision, coagulation, and vaporization of soft and fibrous tissue, including osseous tissue



OtoProbe™

- Available in both long- and short-angled styles
- Bayonet designs for easy viewing of target tissue
- Short-angled OtoProbe provides better view of tissue around the needle tip.
- Long-angled OtoProbe offers a longer, straight needle for easier tissue appplanation



FlexFiber™

- Laser fiber for effective treatment of polyps and adhesions
- Flexible fiber tip for easy access to the target site via a cannula
- Fiber sizes available in 200, 300, 400 and 600 μm diameters





ENT Laser Delivery Devices

Specifications

Compatible Laser System:	OcuLight® TX (2.5 W max) OcuLight GLx (1.5 W max)
Packaging:	6 per box Sterilized, Single Use
Fiber Length:	2.5 m

Product Information

Product Name	Description	Part Number
OtoProbe™	Long-Angled	14310
OtoProbe	Short-Angled	14320
FlexFiber™	200 µm Laser Fiber	15702
FlexFiber	300 µm Laser Fiber	15703
FlexFiber	400 µm Laser Fiber	15704
FlexFiber	600 µm Laser Fiber	15706



Specifications are subject to change without notice. IRIDEX, the IRIDEX logo, and OcuLight are registered trademarks and OtoProbe is a trademark of IRIDEX Corporation.

Products are covered by one or more of the following U.S. patents: 5,085,492; 5,088,803; 5,372,595; 5,511,085; 5,982,789; 6,327,291; 6,540,391; and 6,733,940.

Other U.S. and International patents pending.



IRIDEX