



Eye Safety Filters

Something special
here



iridex.com

Something unique about our filters here

- 532 nm, 577 nm, and 810 nm Eye Safety Filters for Endophotocoagulation
- Superior color balance and clear tissue visualization with 100 coating layers
- Compatible with standard Zeiss, Leica/Wild, Topcon, and Moller Wedel/Haag-Streit microscopes
- Fixed, 2-channel and 4-channel models available
- V notch filter design, low profile housing

Part Numbers and Compatibility	
Zeiss	10378, 30367, 30493, 65471, 65473, 65475
Leica/Wild	11175, 30364, 65470, 30494, 65472, 65474
Topcon	11326, 10378, 30366, 30367, 30493, 65473, 65475, 65551
Moller Wedel/Haag-Streit	13107, 30775, 66076
Alcon	10378, 30367, 65471, 30493, 65473, 65475
Takagi	10378, 30367, 65471, 30493, 65473, 65475
All Iridex 810 nm*	13107, 10378, 11175, 11326, 14080
TX/GLx/GL/IQ532*	30364, 30366, 30367, 30775, 31260
IQ 577*	65470, 65471, 65551, 66076,
TX/GLx/GL**	30493, 30494
IQ 577**	65472, 65473
IQ 532**	30493, 30494, 65474, 65475

* Iridex Fixed Eye Safety Filter.

** Iridex 2-Position Eye Safety Filter Integrated SmartKey Cable connects to console

EC

REP

CE

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Contact Iridex today to learn more
+1.650.962.8100 · customerservice@iridex.com · iridex.com



© 2024 Iridex. All rights reserved.
Iridex and the Iridex logo are registered trademarks of Iridex.
88291 Rev. B 06.2024